

INFORMATION SHEET

REPRODUCTION MOREL MUSHROOMS PAINT SCHEDULE

These artificial mushroom when painted, will be hard to distinguish between the real ones. Just follow this easy paint schedule!

PRE-PAINTING PREPARATION:

To prepare, wipe down the reproduction Morel mushroom with lacquer thinner, paying special attention to the fringes of the mushroom. You can use a variety of colors if you choose. Three colors that can be used and create great results are:

- Yellow Ochre
- Buckskin Tan
- Dark Brown

It is highly recommended that you find some clear reference photos of the mushroom in the stage of growth that you are trying to recreate and study it. Copy the colors and refer to these pictures while painting for a realistic reproduction Morel mushroom!

PAINTING PROCEDURE:

1. Using your airbrush with the color you have chosen, Rotate the mushroom in one hand, carefully misting the edges of the Morel. Work your way completely around one way and then using the same method, airbrush in the opposite direction. The fog from the overspray will work its way down onto the body of the mushroom. You can also mist the stem of the morel mushroom.
3. Now, once the paint you have chosen has dried, use a topcoat such as a base sealer or a gloss color if you choose. The color of the mushrooms will vary depending on the development stage of the specimen that you are trying to recreate.

Have fun and be creative, your sure to have great results!

RESEARCH MANNIKINS

1-800-826-0654 • (541) 451-1538 OFFICE • (541) 451-5455 FAX

“Your Hometown Taxidermy Supply Store”